

SPECIALIZED ACCESSORIES

PerkinElmer's specialized accessories broaden the capability of its spectrophotometers, giving researchers greater control over complex optical measurements. They're designed to deliver consistent, accurate results when characterizing surfaces, coatings, and other challenging materials. By improving precision, stability, and flexibility, these tools make it easier to adapt the instrument to specific experimental needs. The result is reliable, high-quality data that supports deeper insight and more confident interpretation across a wide range of applications.



Absolute Directional VW Specular Reflectance Accessory

L6310205

The absolute directional VW specular reflectance accessory is intended for flat specular samples with a reflectance above 25% such as reflectance standards, laser mirrors, optical solar reflectors and beam splitters. A Glan-Thompson (250 - 2500 nm) polarizer is included. This accessory allows for the calibration of reference mirrors with a reflectance above 80%, absorption measurements of laser mirrors and characterization of beam splitters.